

## 428605

### **QBC Autoread Plus Hematology System**

The QBC Autoread Plus is an affordable and durable dry hematology system that can deliver a 9-parameter Complete Blood Count using QBC Diagnostics' line of unique, acridine orange-coated tubes. The Autoread Plus system comes with the Autoread Plus analyzer, QBC Capillary Centrifuge, Printer (U.S. Only), USB cable, pipetter, AccuTube starter kit, software cartridge, workstation, calibration rod, and power packs for all equipment, including U.S. and European line cords for 120/230 V, 50/60 Hz operation.

## Tubes

### 423406

#### **QBC AccuTubes (Box of 100)**

AccuTubes can be used for capillary samples from a finger or heel stick, as well as venous blood samples. They are packaged with a ventilated stopper inserted into the tube, and a soft grip float holder for easy insertion of the plastic float.

### 424641

#### **QBC E-Z Prep Capillary Tubes (Box of 100)**

E-Z Prep Capillary Tubes are packaged with a ventilated stopper and a pre-inserted float within the tube. Also includes a plastic stylus to ensure proper float placement.

### 424240

#### **QBC Venous Tubes (Box of 100)**

These tubes can be used with the pipetter to draw samples from venous tubes. Contains 10 individual trays of 10 tubes with external floats and manual closures.

### 424237

#### **QBC Venous Tubes (Box of 1000)**

### 424241

#### **QBC Capillary Tubes (Box of 100)**

These tubes can be used with blood samples taken from a finger or heel stick. Contains 10 individual trays of 10 tubes with external floats and manual closures.

### 424238

#### **QBC Capillary Tubes (Box of 1000)**

## Consumables and Accessories

### 424304

#### **Hematology Controls**

These control samples are designed to test abnormal blood cell count values and are a necessity for U.S. customers to meet CLIA requirements for the Autoread Plus. The controls are shipped every three weeks, and contain 6 vials each of 2 different blood samples.

### 424305

#### **Extended Range Controls**

These control samples are designed to meet CLIA requirements for calibration verification. They contain 2 vials each of 3 different blood samples. Note: Unlike hematology controls (424304), extended range controls are available on an as-needed basis.

### 420000

#### **1.8 mm Blade Lancets (Box of 100)**

These spring loaded lancets are ideal for drawing capillary blood for use with QBC Tubes. The green 1.8 mm version is designed for pediatric use.

### 420001

#### **2.3 mm Blade Lancets (Box of 100)**

The blue 2.3 mm Blade Lancets are designed for use on adult patients.

## **Replacement Parts**

**429577**

### **QBC Autoread Plus Analyzer with Accessories**

Autoread Plus analyzer with software cartridge, workstation, AccuTube starter kit, power pack, and U.S. and European line cords for 120/230 V, 50/60 Hz operation. Note: Does not contain centrifuge or printer.

**429576**

### **QBC Autoread Plus Analyzer**

Autoread Plus tube analyzer without accessories, centrifuge, or printer for 120/230 V, 50/60 Hz operation.

**425740**

### **QBC Capillary Centrifuge**

This centrifuge can be used to spin all QBC tubes (except for the STAR Tube) for testing in the Autoread Plus analyzer. The centrifuge has 20 tube slots, and rotates at 12,000 RPM. Comes with power pack and U.S. and European line cords for 120/230 V, 50/60 Hz operation.

**421550**

### **USB Cable**

This 6-foot long cable is required to connect the Autoread Plus analyzer to an external USB-compatible printer.

**424613**

### **Calibration Check Rod**

This specially coded metal carrier rod is used to calibrate the Autoread Plus.

**424226**

### **QBC Workstation**

Sturdy plastic workstation that can be used to hold 100 tubes of varying sizes, including QBC tubes.

**424225**

### **QBC Pipetter**

This pipetter is used with the AccuTube or QBC Venous Tube to draw the correct amount of blood into the tube for testing. Includes AccuTube Pipetter Spacer.

**421705**

### **Pipetter Barrel**

Replacement barrel (the portion of the pipetter into which tubes are inserted) for the pipetter.

**421706**

### **QBC AccuTube Pipetter Spacer**

This spacer is attached to the pipetter to draw the correct amount of blood into the AccuTube during venous tests.

**424593**

### **QBC Autoread Plus Operator's Manual**

A replacement copy of the Autoread Plus manual.

**421923**

### **QBC Autoread Plus Analyzer Software Cartridge**

This replacement cartridge serves as the central processing unit and provides USB printer capacity for the Autoread Plus analyzer. It can be plugged into an Autoread Plus analyzer with a broken or defective software cartridge.

**424590**

### **QBC Autoread Plus Analyzer Power Pack**

Replacement power pack for the Autoread Plus tube analyzer, which can be used with any QBC line cord for 120/230 V, 50/60 Hz operation.

**421763**

### **QBC Capillary Centrifuge Power Pack**

Replacement power pack for the Capillary Centrifuge, to be used with any QBC line cord for 120/230 V, 50/60 Hz operation.

**421634****QBC Line Cord (USA)**

Replacement line cord for the Autoread Plus tube analyzer or Capillary Centrifuge for use in the United States.

**421554****QBC Line Cord (U.K.)****421551****QBC Line Cord (Europe)****421291****QBC Capillary Centrifuge Cover Assembly**

Replacement cover for the inside of the Capillary Centrifuge. Contains both the metal disc cover and the rotor cover access knob.

**421317****QBC Capillary Centrifuge Rotor Cover Access Knob**

Replacement rotor cover access knob for the inside of the Capillary Centrifuge. Contains only the rotor cover access knob.



168 Bradford Drive  
Port Matilda, PA 16870  
814-692-7661  
[www.qbcdiagnostics.com](http://www.qbcdiagnostics.com)