

429000

QBC STAR™ Hematology Analyzer

The QBC STAR is a dry hematology analyzer that delivers a 9-parameter Complete Blood Count with the touch of a button. The STAR's unique microprocessor performs self-calibration, ensuring accuracy. The STAR comes with 2 replacement 3.15 A fuses, 2 rolls of thermal printer paper, screwdriver, forceps, and power cord for 120 V, 60 Hz operation.

Tubes

429625

QBC STAR Tubes (Box of 100)

A box of 100 STAR Tubes. STAR Tubes are the only tubes designed to work with the STAR.

429011

QBC STAR Tubes (Bag of 20)

Consumables & Accessories

424304

Hematology Controls

These control samples are designed to test abnormal blood cell count values, and are a necessity for U.S. customers to meet CLIA requirements for the STAR. The controls are shipped every three weeks, and contain 6 vials each of two different blood samples.

424305

Extended Range Controls

These control samples are designed to meet CLIA requirements for calibration verification. They contain 2 vials each of 3 different blood samples. Note: Unlike hematology controls (424304), extended range controls are available on an as-needed basis.

420000

1.8 mm Blade Lancets (Box of 100)

These spring loaded lancets are ideal for drawing blood for use with the STAR Tubes. The green 1.8 mm version is designed for pediatric use.

420001

2.3 mm Blade Lancets (Box of 100)

The blue 2.3 mm Blade Lancets are designed for use on adult patients.

429580

Printer Paper (3/pk)

Three rolls of 2.5 inch thermal paper for use with the STAR's internal printer.

429581

Fuses, 250 V 3.15 A 5x20 mm (5/pk)

Package of five replacement fuses for the STAR.

429640

Patient ID Upgrade Kit

The kit contains a keyboard that can be used to enter Patient and Operator ID numbers into the STAR.

421634

QBC Line Cord (USA)

Replacement line cord for STAR analyzers in use within the United States.

421554

QBC Line Cord (U.K.)

Replacement line cord for STAR analyzers in use within the United Kingdom.

421551

QBC Line Cord (Europe)

Replacement line cord for STAR analyzers in use within Europe.

429627

Carrying Case

This hard plastic case comes with wheels and foam inserts to easily transport and protect the STAR during travel.



168 Bradford Drive
Port Matilda, PA 16870 U.S.A.
+1-814-692-7661
www.qbcdiagnostics.com